



RCPN

## Standard operating procedure of Nanodrop

**Objective:** This instruction explains how use and keep the Nanodrop.

**Scope:** using for documenting the concentration of nucleic acids or proteins.

**Responsibility:** Laboratory assistance is responsible for the accuracy of the apparatus and should be reported to him/her in case of problems.

**Materials and equipment:** ---

### Guidelines for use:

- 1- Turn on the Nanodrop and its computer.
- 2- Run the ND-1000 software on the desktop.
- 3- Based on your sample type, select Nucleic acid or protein options.
- 4- Wash the sensor of the apparatus with 2  $\mu$ L deionized water (DW) and whatman filter paper.
- 5- Again, add 2  $\mu$ L deionized water on the center of sensor, lower down the arm of Nanodrop and press the “blank” key in the software.
- 6- Again remove the DW and add 2  $\mu$ L of your source solution that used for extraction procedure on the sensor and press the Blank key.
- 7- After removing the source solution, add 2  $\mu$ L of your sample and lower down the arm and press the “Measure” key.
- 8- Document the results and after washing the sensor with DW, exit from the software and turn off the apparatus.

### Warning and safety precautions:

If any abnormal sound or vibrate are felt, press the “Stop” bottom immediately.

### Maintenance:

The door of the apparatus should be always closed.

Document ID: RCPN 46

Document name	Standard operating procedure of Nanodrop
Document arrangement	Mohanna Osali
Date	10/25/2018
Full file name	Standard operating procedure of Nanodrop
Document Description	This instruction explains how use and keep the Nanodrop.
Author / translator	Nazli Saeedi